

Instructions for Use

Prior to Use:

1. Perform perineal skin assessment and document per hospital protocol.
2. Obtain necessary supplies: The CAREDRY® Male device, one or two lengths of tubing (depending on mounted or free-standing canister), and a suction canister.

Setup:

1. Using a length of tubing, connect a clean collection canister to a suction regulator that is set to 125 mmHg.
2. Attach an additional length of tubing to the canister, turn suction to CONTINUOUS and listen for airflow.
3. Remove the CAREDRY® Male device from its packaging and connect to tubing and collection canister.
4. Peel off the back liner using the pull tab to expose the adhesive pad for placement.
5. With the narrow end of the bag pointed down (towards the patient's legs), place the patient's penis through the hole in the adhesive pad so that the scrotum remains external, and that the penis is fully contained within the pouch. Once in place, gentle pressure should be applied around the opening to secure the device.
6. Ensure that the penis can be visualized through the clear front panel of the device and that it is generally pointed towards the narrow end of the device - where suction is applied.

In Use:

1. Continue to assess skin condition and ensure tubing is connected while in use.
2. Assess device periodically to ensure tubing is not a discomfort to patient
3. Replace when soiled, and every two days or per hospital protocol.
4. Replace suction tubing and empty and/or replace collection canister as needed per hospital protocol.

Replacement:

1. Obtain replacement device.
2. Turn suction to OFF.
3. Gently lift the adhesive pad to remove device and dispose of the unit.*
4. Assess skin integrity & document per hospital protocol.
5. Using a new System, repeat the setup instructions, beginning at step 3.

* System may be considered a potential biohazard. Dispose of according to hospital protocols.

BOEHRINGER

®

Instructions for Use

CAREDRY® MALE

Description

CAREDRY® Male is a non-invasive urinary management system. The System features a flexible pouch, gentle adhesive and propriety venting system.

Indications

CAREDRY® Male is indicated for external urinary incontinence management.

Contraindications

Do not use for applications other than collection of urine, or on patients with normal urinary function.

Do not use on patients with bowel incontinence without adequate management to reduce risk of contamination.

Disposal

Dispose per hospital protocol or sooner if there is evidence of contamination.

Empty or replace collection canister as needed, per hospital protocol.

CAREDRY® Male should be disposed of per your facility's bio-hazardous waste protocols.

Storage & Shelf Life

These devices are manufactured with a 1 year shelf life. Always store in a climate controlled space between 60 and 80 degrees Fahrenheit.

Warnings



CAREDRY® Male is for single patient use. Reuse may adversely affect the performance and risks cross contamination of patients. Do not attempt to clean.



CAREDRY® Male is packaged clean; but non-sterile



Latex Free



CAREDRY® System must only be attached to vacuum systems. Do not attach to compressed air, nitrogen, or oxygen.

Ensure wall outlet meets the requirements of NFPA 99. Non-compliant wall outlets may affect regulator performance.

Discontinue use if skin breakdown or allergic reaction occurs.

Do not insert into any body cavity.

Do not use with portable suction pumps.

Do not attach to intermittent suction.

Do not use with unregulated suction.

Do not use with bedpan or other materials that cause insufficient airflow to the device.



Boehringer Laboratories, LLC.
300 Thoms Drive, Phoenixville, PA 19460
1-800-642-4945 www.boehringerlabs.com
Made in USA

BOEHRINGER

®

Caring for Lives through Innovation, Quality and Service

Part Number	Description	Ordering
7680	CAREDRY® Male	1 Box (20 Sytems)